

SARS-CoV-2 Surveillance Shipping Manifest



Facility:

Phone Number:

Sender Name:

Address:

Date:

Samples can be refrigerated (2-8 °C) and shipped with cold packs if received by the PHL within 72 hours of collection. Any specimens that will not arrive at the PHL within 72 hours of collection must be stored frozen (-70 °C or colder) and shipped on dry ice, along with this completed manifest. You may batch samples and submit them to the PHL once per week.

Please list out all samples in this shipment:

Collection Date	Patient Name	Sex	Date of Birth	City, State	Zip